



# Planmeca Romexis® 3D Ortho Studio



# 3D tools for orthodontists and dental labs

All digital dental models in STL format can be used in Planmeca Romexis® 3D Ortho Studio

The **Planmeca Romexis® 3D Ortho Studio** module offers innovative tools for orthodontists and dental laboratories. It is designed for the examination and analysis of scanned digital dental models, as well as for orthodontic treatment planning in 3D.



## Alignment of digital dental models

Digitise dental models by scanning dental impressions and plaster casts with the **Planmeca ProMax® 3D** model scanning program. Alternatively, direct digital impressions can be taken with the **Planmeca Emerald™** or **Planmeca PlanScan®** intraoral scanners. Examination, analysis and treatment planning can be done conveniently in the **Planmeca Romexis® 3D Ortho Studio** module.



Plaster cast in Planmeca ProMax® 3D



Impression scan in Planmeca ProMax® 3D




Planmeca Emerald™ full arch scan

## Preparation of digital study models

Utilise the module's efficient sculpting and smoothing tools for removing artifacts and excess material from digital dental models. A virtual base can be attached to the digital models to make them more presentable to the patient and more suitable for analysis.

## Dental model analyses

Planmeca Romexis 3D Ortho Studio includes versatile tools for dental model analysis. Tooth, arch and cross-sectional measurements can be used for tooth size, space, overjet and overbite analyses. In addition, teeth intersections and occlusions can be visualised with a deviation colour map.



## Treatment planning in the virtual setup

In the virtual setup, a staged treatment plan can be established by displacing teeth while visualising intersections and contacts. Attachments and brackets can be placed on top of the teeth to create indirect bonding templates.

Animations of tooth movement visualise the treatment progression from initial setup to the planned treatment outcome, and improves communication between the orthodontist and patient.



Digital model → Digital aligner → 3D printed aligner → 3D printed dental arch → Thermoformed aligner

## Creating and exporting a digital model series

**Planmeca Romexis® 3D Ortho Studio** automatically generates an incrementally progressing series of digital models from the initial setup to the planned treatment outcome based on user given parameters for maximum tooth movement and rotation between models. Virtual clear aligner models can be automatically created on top of the digital models.

All the models can be exported in STL format for 3D printing, custom appliance design and manufacturing. Excess material, such as the virtual base, can be cut off from the models before export.

## Extensive analysis and treatment plan report

All the measurements and treatment planning information, such as tooth movements, stripping and extractions, can be summarised in a report together with model series information.

## Planmeca Romexis® 3D Ortho Studio Viewer for sharing data with partners

Share an ortho project with a patient with the free **Planmeca Romexis® 3D Ortho Studio Viewer**, for example by using the **Planmeca Romexis® Cloud** image transfer service.

Exported digital models can be printed with any compatible 3D printer or manufacturing machine

The module is compatible with the Windows operating system.

# PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentreVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca CALM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca Creo, Planmeca Emerald, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Olo, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanID, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca SmartGUI, Planmeca Solanna, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Viso, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy, and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

