

# Planmeca Compact™ i Classic

The ideal choice for the needs of general dentistry



Simply  
efficient

# Optimal access with a good-looking design

The functional design of **Planmecca Compact™ i Classic** guarantees a smooth workflow and an excellent working position.

## Timeless design, robust construction

We use aluminium parts and an epoxy paint finish in our dental units for ultimate durability. The seamless, powder-coated surfaces meet the strictest infection control requirements.

## Perfect access

The side-lifting chair and small cuspidor base maximise the legroom for both dentist and assistant. The lean and narrow backrest design provides excellent access to the treatment area.

## Unique linkage of dental unit and software

**Planmecca Romexis® Clinic Management** can connect all your dental units to a network to generate valuable data of dental unit usage and events for efficiency and quality assurance purposes.



# Easy to learn, easy to use

*Planmecca Compact™ i Classic has a traditional keyboard with all the essential unit, instrument and infection control functions. A keyboard is also available for the assistant suction arm for easy and ergonomic use.*

## Perfectly balanced instrument arms

The balanced instrument arms guarantee ergonomic instrument use since only a small movement is sufficient to activate the instrument. Alternatively, **Planmecca Compact™ i Classic** can also be equipped with hanging tube instruments.

## Instant adaptation to different needs

The freely interchangeable instruments and easy console movements enable swift conversion to personal preferences without any tools.



## Intraoral scanner just like any other instrument

The seamless dental unit integration of the intraoral scanner makes real-time scanning data immediately available from your chairside tablet, guaranteeing a smooth workflow at all times.

## All operations from one foot control

All the necessary functions can be carried out from the unit's robust foot control, keeping your hands free and minimising cross-contamination.

# Excellent patient comfort

*The dental chair design follows the shape of the patient's body and provides firm and comfortable support.*

## Ergonomically designed chair

The double articulated headrest allows comfortable positioning of the patient's head. Thanks to the smooth headrest length adjustment, the chair adapts easily to patients of all sizes.

Our dental units can also be fitted with armrests for additional patient support.

## High-quality upholstery for enduring comfort

The **Comfy™** upholstery is made of durable artificial leather. The thin shape ensures excellent reach to the treatment area while the seamless design allows easy upkeep and enhanced infection control.

Made of extremely durable artificial leather and viscoelastic memory foam, the **Ultra Relax™** upholstery adapts perfectly to the patients' body shape, allowing them to stay relaxed even during the longest treatment sessions.



Child cushion set offers optimal support for smaller patients

<b>Fuchsia</b> 6583 • Ultra Relax 710 • Comfy	<b>Cassis</b> 9314 • Ultra Relax	<b>Purple</b> 9333 • Ultra Relax 706 • Comfy	<b>Red</b> 1176 • Ultra Relax 702 • Comfy
<b>Metallic Silver</b> 5002 • Ultra Relax	<b>Sun</b> 5250 • Ultra Relax	<b>Lime</b> 4482 • Ultra Relax 602 • Comfy	<b>Petrol</b> 2520 • Ultra Relax 701 • Comfy
<b>Sky</b> 2665 • Ultra Relax 707 • Comfy	<b>Azure</b> 2674 • Ultra Relax 708 • Comfy	<b>Deep Sea Blue</b> 604 • Comfy	<b>Metallic Blue</b> 4020 • Ultra Relax
<b>Steel Blue</b> 401 • Comfy	<b>Baltic</b> 2551 • Ultra Relax	<b>Dark Blue</b> 504 • Comfy	<b>Diplomat Blue</b> 2478 • Ultra Relax
<b>Metallic Gold</b> 3174 • Ultra Relax	<b>Hazel</b> 704 • Comfy	<b>Taupe</b> 3779 • Ultra Relax	<b>Walnut</b> 3777 • Ultra Relax
<b>Coffee</b> 705 • Comfy	<b>Graphite</b> 5806 • Ultra Relax 703 • Comfy	<b>Black</b> 5607 • Ultra Relax 305 • Comfy	



One software for all.



*Planmeca Oy designs and manufactures a full line of industry-leading dental equipment, including 3D and 2D imaging devices, CAD/CAM solutions, dental care units and software. Planmeca Oy, the parent company of the Finnish Planmeca Group, is strongly committed to better care through innovation, and it is the largest privately held company in the field.*

# PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentreVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca CALM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca Creo, Planmeca Emerald, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Olo, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanID, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca SmartGUI, Planmeca Solanna, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Viso, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy, and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

10030886/1017/en

