

EXPLORERS

For diagnosing calculus and caries and exploring of pockets, restorations and furcations.



Explorer 5-8

- LMErgoSense LM 5-8 ES
- LMErgoMax LM 5-8 XSI
- LMErgoNorm LM 5-8 SI

- Standard explorer for diagnostic and operative procedures



Explorer 5-9

- LMErgoSense LM 5-9 ES
- LMErgoMax LM 5-9 XSI
- LMErgoNorm LM 5-9 SI

- Standard explorer for diagnostic and operative procedures



Especially for probing distal surfaces of posterior teeth.

Explorer, centered 9-17C

- LMErgoSense LM 9-17C ES
- LMErgoMax LM 9-17C XSI
- LMErgoNorm LM 9-17C SI

- Centered tip (17C)
- Elongated lower shank and "back-action" blade



Explorer 17-23

- LMErgoSense LM 17-23 ES
- LMErgoMax LM 17-23 XSI
- LMErgoNorm LM 17-23 SI

- For distal surfaces (17)
- Elongated lower shank and "back-action" blade (17)
- For areas in the vicinity of proximal contact surfaces (17)
- Also for probing occlusal surface grooves (23)
- For detecting subgingival calculus (23)



Especially for probing distal surfaces.

Explorer 17C-23

- LMErgoSense LM 17C-23 ES
- LMErgoMax LM 17C-23 XSI
- LMErgoNorm LM 17C-23 SI

- For distal surfaces (17C)
- Centered tip (17C)
- For areas in the vicinity of proximal contact surfaces (17C)
- Also for probing occlusal surface grooves (23)
- For detecting subgingival calculus (23)



Explorer 3-6

- LMErgoSense LM 3-6 ES
- LMErgoMax LM 3-6 XSI
- LMErgoNorm LM 3-6 SI

- Standard explorer for diagnostic and operative procedures



For subgingival root examination.

Explorer 11-12

- LMErgoSense LM 11-12 ES
- LMErgoMax LM 11-12 XSI
- LMErgoNorm LM 11-12 SI

- Angles similar to Gracey 11/12
- For probing proximal and cervical caries



High sensitivity at seams.

Handle options:



LMErgoSense



LMErgoMax



LMErgoNorm

Flexplorer 11F-12F

LMErgoSense LM 11F-12F ES

LMErgoMax LM 11F-12F XSI

- Flexible explorer angulated similarly to Gracey 11/12
- Well-suited for periodontal diagnosing of subgingival calculus and probing deep pockets

Explorer 18-19

LMErgoSense LM 18-19 ES

LMErgoMax LM 18-19 XSI

LMErgoNorm LM 18-19 SI

- For prosthetic crowns, bridges and inlay- / onlay fillings
- Excellent sensitivity at seams

Explorer 29

LMErgoSense LM 29 ES

LMErgoMax LM 29 XSI

LMErgoNorm LM 29 SI

- Holst design
- Lateral angle on the shank
- Single-ended

Explorer 9

LMErgoSense LM 9 ES

LMErgoMax LM 9 XSI

LMErgoNorm LM 9 SI

- Single-ended