

MIRROR HANDLES

Simple stem

LM 25 mirror handle threads are in accordance with international standard M2.5. The mirror shank must have a diameter of 2.5 mm.

Mirrors are sold separately. Product no 10P12, box of 12 pcs, mirror #4 simple stem.

Cone socket

LM 28 mirror handle threads are in accordance with international standard UNF4-48. The mirror shank must have a diameter of 2.6 mm.



ErgoSingle Mirror Handle

- LMErgoSingle** LM 25 ESI
- LMErgoSingle** LM 28 ESI

- Well balanced, better reachability
- Colors: red, grey, blue



Mirror Handle 25

- LMErgoSense** LM 25 ES
- LMErgoMax** LM 25 XSI
- LMErgoNorm** LM 25 SI

- Single-ended
- Simple stem design
- Colors: red, grey, blue



Mirror Handle 28

- LMErgoSense** LM 28 ES
- LMErgoMax** LM 28 XSI

- Single-ended
- Cone socket design, UNF4-48, with special threads
- Colors: red, grey, blue



Mirror handle, ortho

- LMErgoSense** LM 25-26 ES
- LMErgoMax** LM 25-26 XSI
- LMErgoSense** LM 28-26 ES
- LMErgoMax** LM 28-26 XSI

- Bausch-Verbiest design
- Orthodontic mirror handle
- For measuring the sagittal overjet and vertical overbite
- Millimeter scale, total length 25 mm



Mirror Handle JM

- LMErgoSense** LM 25-473 ES
- LMErgoMax** LM 25-473 XSI

- For implant installation

Handle options:

- LMErgoSense**
- LMErgoMax**
- LMErgoNorm**