

LM DarkDiamond

LM

feel the
difference



Pleasant and precise sculpting of composite



Diamond like feel & precision

Restorative instruments with LM Dark Diamond non-stick coating are ideal for placing and contouring of composite resins. The diamond like coating reduces the stickiness and allows pleasant and precise restorative instrumentation.

- Exclusive diamond like coating for non-stick placement and sculpting of composite resin
- Not sensitive to scratches due to hard surface
- Smooth surface for easy cleaning
- Dark color coating offering maximum contrast against filling and minimal reflections
- Ergonomic LM-ErgoSense for excellent precision and feel



Scan the QRs for video and more information

Instrument range



Burnisher

LMErgoSense

LM 31-32DD ES

- LM Dark Diamond non-stick coating
- For packing and modeling composite



Burnisher-Plugger,
ball-ended 1.5-2.5 mm

LMErgoSense

LM 471-473DD ES

- LM Dark Diamond non-stick coating
- For packing and modeling composite



Modeller MH

LMErgoSense

LM 708-709DD ES

- LM Dark Diamond non-stick coating
- For packing large approximal fillings
- Cylinder shaped working end will not get stuck into the composite while light curing



Modeller MH Mini

LMErgoSense

LM 707-708DD ES

- LM Dark Diamond non-stick coating
- For packing small and medium approximal fillings
- Cylinder shaped working end will not get stuck into the composite while light curing



Packer-Modeller IO

LMErgoSense

LM 482-483DD ES

- LM Dark Diamond non-stick coating
- For packing and modeling composites
- Conical working ends optimal for shaping occlusal surfaces



Plastic instrument, anterior

LMErgoSense

LM 43-47DD ES

- LM Dark Diamond non-stick coating
- For placing, packing and contouring composite



Plastic instrument, posterior

LMErgoSense

LM 48-702DD ES

- LM Dark Diamond non-stick coating
- For placing, packing and contouring composite



Mini Spatula

LMErgoSense

LM 444-445DD ES

- LM Dark Diamond non-stick coating
- For placing and contouring composite

LM Dark Diamond instruments can be cleaned and sterilized like any other LM instrument.

LM Dental

LM-Instruments Oy

Norrbyn rantatie 8 | FI-21600 Parainen | Finland

Tel. +358 2 4546 400

info@lm-dental.com | www.lm-dental.com

Contact your local distributor!